2X M5 Taq HiFi PCR mix 使用说明书

产品名称	单位	货号
2×M5 Taq HiFi PCR Mix	1 ml	MF002-01
2×M5 Taq HiFi PCR Mix	5×1 ml	MF002-05
2×M5 Taq HiFi PCR Mix	10×1 ml	MF002-10

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2X M5 Taq HiFi PCR mix

Introduction

2× M5 Taq HiFi PCR mix is a convenient premixed 2× concentrated solution for PCR which includes M5 Taq HiFi DNA Polymerase, PCR Buffer, Mg²⁺, dNTPs, PCR Stabilizer and PCR Intensifier. Utilizing M5 Taq HiFi DNA polymerase with high amplification efficiency and low mismatching rate, the master mix offers the whole reaction system with high sensibility, specificity and stability. It is easy to use and can maximal avoid the personal error and contamination. It is suitable for PCR amplification with high specificity or large scale of gene test of difficult templates with secondary structure. The PCR product has adenine residue at the 3'-end of both strands, and it could be used for TA cloning (Mei5 Bio TOPO-TA cloning kit MF019 or MF020)、enzyme digestion、ligation and sequencing after purification. An inert blue dye (Bromophenol Blue, with migration equivalent to 500bp DNA fragment) is included in the master mix, which allows for direct loading of the PCR product(s) onto an agarose gel.

Components:

Cat. No.	Description	Quantity
MF002-01	2× M5 Taq HiFi PCR Mix	1 ml
MF002-05	2× M5 Taq HiFi PCR Mix	5×1 ml
MF002-10	2× M5 Taq HiFi PCR Mix	10×1 ml

Note: 2×M5 Taq HiFi PCR Mix contains Taq HiFi DNA Polymerase, 300 μM dNTP each, 20 mM Tris-HCl (pH8.7), 100 mM KCl, 3 mM MgCl₂.

Shipping and Storage

Keep at -20°C for long-term storage. 2× M5 Taq HiFi PCR mix are stable at 4°C for six months. For daily use, we recommend keeping an aliquot at 4°C.

Applications:

- 1) Gene test, half-QPCR and detection of trace DNA.
- 2) Gene cloning, SNP analysis and site-directed mutagenesis.

Protocol

PCR reaction system

Template DNA	<1μg
2× M5 Taq HiFi PCR Mix	10 μl
Primer 1 (10 μM)	0.5 μl
Primer 2 (10 μM)	0.5 μl
Nuclease-free ddH ₂ O	Up to 20 μl

PCR reaction conditions

Procedure	Temperature	Time	
Pre-denaturation	94°C	3 min.	} 32-36cycles
Denaturation	94°C	25 sec.	
Annealing	55-64°C	25 sec.	
Extension	72°C	30-60sec. /1 kb DNA	
Final extension	72°C	5 min.	
Keeping	4°C	forever	