## 2× M5 Taq PCR Mix 使用说明书

产品名称	单位	货号
2× M5 Taq PCR Mix	1 ml	MF001-01
2× M5 Taq PCR Mix	5×1 ml	MF001-05
2× M5 Taq PCR Mix	10×1 ml	MF001-10
2× M5 Taq PCR Mix	100×1 ml	MF001-100

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## 2× M5 Taq PCR Mix

### Introduction

Mei5 Bio 2× M5 Taq PCR mix is a convenient premixed 2× concentrated solution for PCR reactions, which includes M5 Taq DNA Polymerase, PCR Buffer, Mg<sup>2+</sup>, dNTPs, PCR Stabilizer and PCR Intensifier. The unique composition offers the whole reaction system with high sensibility, specificity and stability. It is easy to use and can maximal avoid the personal error and contamination. It is suitable for PCR amplification with high specificity or large scale of gene test of complicate templates with secondary structure or GC-rich region. The PCR product has adenine residue at the 3'-end of both strands, and it could be used for TA cloning (Mei5 Bio TOPO-TA cloning kit MF019 or MF020)、enzyme digestion、ligation and sequencing after purification. A red dye (with migration equivalent to 800bp dsDNA) is included in the master mix, which allows for direct loading of the PCR product(s) onto an agarose gel.

### Components

Cat. No.	Description	Quantity
MF001-01	2× M5 Taq PCR Mix	1 ml
MF001-05	2× M5 Taq PCR Mix	5×1 ml
MF001-10	2× M5 Taq PCR Mix	10×1 ml
MF001-100	2× M5 Taq PCR Mix	100×1 ml

**Note:** 2×M5 Taq PCR Mix contains Taq DNA Polymerase, 300 μM dNTP each, 20 mM Tris-HCl (pH8.7), 100 mM KCl, 3 mM MgCl<sub>2</sub>.

### Shipping and Storage

Keep at -20°C for long-term storage. 2× M5 Taq PCR mix are stable at 4°C for six months. For daily use, we recommend keeping an aliquot at 4°C.

### Applications

- 1) Gene test, half- QPCR and detection of trace DNA.
- 2) Gene cloning, SNP analysis and site-directed mutagenesis.

### Protocol

#### PCR reaction system

Template DNA	<1μg
2× M5 Taq PCR Mix	10 μl
Primer 1 (10 μM)	0.5 μl
Primer 2 (10 μM)	0.5 μl
Nuclease-free ddH <sub>2</sub> O	Up to 20 μl

#### PCR reaction conditions

Procedure	Temperature	Time	
Pre-denaturation	94°C	3 min.	
Denaturation	94°C	25 sec.	} 32-36cycles
Annealing	55–64°C	25 sec.	
Extension	72°C	30 sec.-1 min./1 kb DNA	
Final extension	72°C	5 min.	
Keeping	4°C	forever	