2x Realtime PCR Super mix (SYBRgreen,with anti-Taq) 使用说明书

产品名称	单位	货号
2xRealtime PCR Super mix	1ml*5 支	MF013-01
2xRealtime PCR Super mix	(1ml*5 支)*5	MF013-05
2xRealtime PCR Super mix	(1ml*5 支)*10	MF013-10

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2x Realtime PCR Super mix (SYBRgreen, with anti-Taq) User Manual

Product	Units	Cat.#
2xRealtime PCR Super mix	1ml*5tubs	MF013-01
2xRealtime PCR Super mix	(1ml*5tubes)*5	MF013-05
2xRealtime PCR Super mix	(1ml*5tubes)*10	MF013-10

Store at -20°C, protected from light.

Description

MF013 is a Taq DNA polymerase-based 2x Super mix for real-time PCR, which contains all components, except for the primer. This reagent is applicable for intercalation assay with SYBR® Green I. This reagent can be used in glass capillary systems (e.g., LightCycler, Roche Molecular Systems, Inc.). This reagent can be used in a passive reference system (e.g., ABI PRISM® 7700, Applied Biosystems, Inc.). The passive reference dye does not affect any other systems. Hot Start technology with anti-Taq DNA polymerase antibodies enables high specificity and reproducible amplification.

Detection

MF013 can be used in general detection devices, such as: LineGene (Bioer Technology co., ltd.); MF013 can also be used in detection equipment using glass capillaries or passive reference, such as: LightCycler (Roche Molecular Systems); ABI PRISM® 7000, 7700, and 7900 (Applied Biosystems). Note: The passive reference mode of detectors should be set at "ROX".

Specimen

1. cDNA: Reverse transcription reactions from total or poly (A)+ RNA may be used directly, or after dilution, for real-time PCR. Purified cDNA by phenol/chloroform extraction and ethanol precipitation may also be used. Oligo dT and random primers are suitable for the reverse transcription reaction. Up to 20% of the synthesized cDNA solution from the M5 Super qPCR RT kit (Code NO. MF012) may be added to the PCR reaction solution directly, without purification.
2. Genomic DNA: Purified DNA, which would be used for general PCR, is also suitable for real-time PCR. In the case of mammalian genomic DNA, 1~10 ng genomic DNA is sufficient for real-time PCR.

Protocol I
1. Intercalation assay protocol using ABI PRISM® 7700

The following is an intercalator assay protocol to be used with ABI PRISM® 7700. For other detection devices, this protocol may require modification depending on each instruction manual.

(1) Preparation of reaction solution

Component	Volume	Final Concentration
PCR grade water	16 μ l	
Realtime PCR Super Mix	25 μ l	1x
10pmol/ μ l (10 μ M) Primer #1	2 μ l	0.4 μ M
10pmol/ μ l (10 μ M) Primer #2	2 μ l	0.4 μ M
Template DNA	5 μ l	
Total volume	50 μ l	

Notes: The primer concentration can be further optimized, if needed. The optimal range for the primers is 0.2~0.6 μ M. In the case of commercially available primers, recommended conditions from those companies should be used.

(2) Cycling conditions <3-step cycle>

The following condition is recommended:

Pre-denaturation:	95°C, 30sec. ~1 min.	
Denaturation:	95°C, 15 sec.	
Annealing:	55~65°C, 15 sec.	
Extension:	72°C, 45 sec. (data collection)	(35~40 cycles)
Melting curve analysis		

Notes:

- The annealing temperature in 3-step cycle should be set to 55~65°C, depending of the primer T_m value.
- The pre-denaturation condition described above is sufficient for inactivation of the anti-Taq DNA polymerase antibodies used in Hot Start PCR. To prevent unexpected and inappropriate results, do not prolong the pre-denaturation period. Fifteen seconds is also sufficient for denaturation during each cycle.
- Data collection step should be longer than 30 sec.

Protocol II

2. Intercalation assay protocol using Roche LightCycler™

The following is an intercalator assay protocol to be used with the Roche LightCycler™. In the case of other detection devices, this protocol should be modified accordingly.

(1) Preparation of reaction solution

Component	Volume	Final Concentration
PCR grade water	6.4 μ l	
Realtime PCR Super Mix	10 μ l	1x
10pmol/ μ l (10 μ M) Primer #1	0.8 μ l	0.4 μ M
10pmol/ μ l (10 μ M) Primer #2	0.8 μ l	0.4 μ M
Template DNA	2 μ l	
Total volume	20 μ l	

Notes: The primer concentration can be further optimized, if needed. The optimal range for primers is 0.2~0.6 μ M. In the case of commercially available primers, recommended conditions from each manual should be followed.

(2) Cycling conditions <3-step cycle>

The following condition is recommended:

Pre-denaturation:	95°C, 30 sec. ~1 min.	
Denaturation:	95°C, 5 sec.	
Annealing:	55~65°C, 10 sec.	
Extension:	72°C, 15 sec. (data collection)	(35~40 cycles)

Melting curve analysis

Notes:

- The annealing temperature can be set to 55~65°C, depending on the primer T_m value.
- The annealing time should be set for 5~20 seconds. Longer annealing time results in increased efficiency, and a shorter time decreases non-specific amplification.
- The pre-denaturation condition described above is sufficient for inactivation of the anti-Taq DNA polymerase antibodies used in Hot Start PCR. To prevent unexpected and inappropriate results, do not prolong the pre-denaturation period. Five seconds is also sufficient for denaturation during each cycle.
- Data collection step should be longer than 10 sec. If commercially available primers or probes are employed, the recommended conditions from each company should be used.