

## M5 HiPer Freund's Adjuvant, Complete

### 弗氏完全佐剂使用说明书

产品名称	单位	货号
M5 HiPer Freund's Adjuvant, Complete	10ml	MF815-01
M5 HiPer Freund's Adjuvant, Complete	10x10ml	MF815-10

#### 【储存条件】

4°C保存。

#### 【产品简介】

Freund's Adjuvant is one of the most commonly used adjuvants in research. It is used as a water-in-oil emulsion. It is prepared from non-metabolizable oils (paraffin oil and mannide monooleate). If it also contains killed Mycobacterium tuberculosis it is known as Complete Freund's Adjuvant. Without the bacteria it is Incomplete Freund's Adjuvant. First developed by Jules Freund in the 1940's, Freund's Adjuvant is designed to provide continuous release of antigens necessary for stimulating a strong, persistent immune response. The main disadvantage of Freund's Adjuvant is that it can cause granulomas, inflammation at the inoculation site and lesions. The mycobacteria in Complete Freund's attracts macrophages and other cells to the injection site which enhances the immune response. For this reason, the Complete Freund's Adjuvant is used for the initial injections. To minimize side-effects, Incomplete Freund's Adjuvant is used for the boosts.

#### 【产品组份】

Each ml contains 1 mg of heat-killed and dried Mycobacterium tuberculosis (strain H37Ra, ATTC25177), 0.85 ml paraffin oil and 0.15 ml of mannide monooleate.

#### 【操作步骤】

1. Using Complete Freund's Adjuvant, vortex or shake to resuspend the Mycobacterium.
2. Mix antigens (preferably in saline) with an equal volume of the adjuvant to form an emulsion. In order to do this, vigorous and prolonged mixing is needed. There are at least three methods which can be used to accomplish this: For small volumes the emulsion can be made in a tube. Pipet the adjuvant in the tube first. Then, while vortexing, add an equal volume of the antigen solution. Vortex vigorously until a thick emulsion forms.

For intermediate volumes, use two syringes connected through a luer fitting. Ideally, a 3-way valve should be used. Take the desired amount of antigen solution into a glass syringe. The volume should not fill more than half the syringe. Take an equal volume of the adjuvant into another glass syringe. Remove all air and connect the syringes through the luer fitting to the 3-way valve. Adjust the 3-way valve such that the connection is open between the two syringes. Carefully depress the plunger from the antigen solution first, pushing the antigen into the oil of the adjuvant. Alternately push the plungers, mixing the adjuvant and the antigen solution into an emulsion. Continue until the plungers are difficult to push.

For large volumes, use a tissue homogenizer to make the emulsion. Add the adjuvant to the homogenizer first. Run the homogenizer for a short time to coat the inside with the adjuvant. Add an equal volume of the antigen solution and run until a thick emulsion forms.

3. The resulting emulsion should be very thick and a drop of it should not disperse if tested by placing on the surface of a saline solution.
4. Transfer the emulsion to a syringe (or remove one syringe from the luer fitting if using the two-syringe method). Remove all the air. Add an appropriately sized needle. The samples are now ready for injection.

**【注意事项】**

These products are for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

**【备注】**

本产品仅供科研使用。在确认产品质量出现问题时，本公司承诺为客户免费更换等量的质量合格产品。