# M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa) 超宽范围 Western 曝光蛋白 Marker 使用说明书

产品名称	单位	货号
M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa)	250 $\mu$ l	MF291-01
M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa)	2x250 $\mu$ l	MF291-02

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M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa) is a ready-to-use mixture with 10 IgG-binding proteins covering a wide range of molecular weights from 15 to 200 kDa in Tris Glycine buffer.

M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa) performs dual functions. First, it contains 4 pre-stained proteins (10, 25, 45 and 70 kDa) for monitoring protein separation during SDS-PAGE, verification of Western transfer efficiency on membranes (nitrocellulose, PVDF, or nylon) and for approximating the protein size. Second, ten IgG binding proteins can be immuno-detected on film or by CCD imaging.

M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa) is compatible for chemiluminescent, fluorescent, chromogenic or other detection systems. In addition, M5 SuperRange Self-Lighting Western Protein Marker (15-200 kDa) has two reference bands with enhanced intensity (at 30 kDa and 80 kDa). The marker is supplied in the gel loading buffer and is ready to use. Do NOT heat, dilute, or add reducing agents before loading.

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
- a. Transfer with buffer containing 20% methanol to fix proteins on membrane.
- b. Wash membrane with PBS or TBS containing less than 0.1% Tween-20 at 4°C.

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1X Tris buffered saline: 6.05 g Tris and 8.76 g NaCl in 800 mL of H<sub>2</sub>O.  
Adjust pH to 7.5 with 1 M HCl and make volume up to 1 L with H<sub>2</sub>O

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