M5 Prestained Plus Protein Ladder (10-250 kDa) with 40kD 使用说明书

产品名称	单位	货号
M5 Prestained Plus Protein Ladder with 40kD	250 μ l	MF028-plus-01
M5 Prestained Plus Protein Ladder with 40kD	2x250 μ l	MF028-plus-02

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M5 Prestained Plus Protein Ladder (10-250 kDa) with 40 kD User Manual

Product Name	Units	Cat.#
M5 Prestained Plus Protein Ladder with 40kD	250μl	MF028-plus-01
M5 Prestained Plus Protein Ladder with 40kD	2x250μl	MF028-plus-02

Storage: Upon receipt store at -20°C. Product is shipped with an ice pack.

Storage Buffer: 62.5mM Tris•H₃PO₄ (pH 7.5 at 25°C), 1mM EDTA, 2% (w/v) SDS, 10mM DTT, 1mM NaN₃ and 33% (v/v) glycerol.

Introduction:


The M5 Prestained Plus Protein Ladder is a prestained mixture of eleven recombinant proteins ranging from 10kDa to 250kDa. Two different chromophores are bound to the proteins, producing a brightly colored ladder (see website for product images). The protein ladder is conveniently packaged and ready to use with no heating, diluting or additional reducing agent necessary.

Important Product Information

- Do not boil the protein ladder.
- The molecular weights of the proteins have a lot-to-lot variation of approximately 3%.
- In low-percentage gels (< 10%), the low-molecular weight proteins in the ladder may migrate with the dye front.
- The large proteins (> 100K) in the ladder may require longer transfer times or higher transfer voltages for Western blotting.
- The mobility of prestained proteins can vary in different SDS-PAGE buffer systems; however, they are suitable for approximate molecular weight determination when calibrated against unstained standards in the same system.

Procedure for Use in Polyacrylamide Gel Electrophoresis:

1. Thaw the ladder at room temperature. Do not boil protein ladder.
2. Mix gently and thoroughly to ensure that the solution is homogeneous.
3. Load an appropriate volume of the ladder onto the gel.
 - Mini-gel: 5μL per well (0.75-1.0mm thick) or 10μL per well (1.5mm thick)
 - Large gel: 10μL per well (0.75-1.0mm thick) or 20μL per well (1.5mm thick)
4. Return the unused protein ladder to -20°C for up to one year or 4°C for up to three months.

kDa, 4-15% PAGE