

## M5 Rabbit Anti-Goat IgG, FITC Conjugated

Product	Unit	Cat.#
M5 Rabbit Anti-Goat IgG, FITC Conjugated	100ul	MF262-01

**【Storage】** : Store at **-20°C** for one year. Avoid repeated freeze/thaw cycles.

英文名称 Rabbit Anti-Goat IgG H&L/FITC

别 名 Rabbit Anti-Goat IgG H&L (FITC); Immunoglobulin G;

抗体来源 Rabbit

克隆类型 Polyclonal

交叉反应 Goat,

产品应用 Flow-Cyt=1:100-1000 IF=1:100-1000

not yet tested in other applications.  
optimal dilutions/concentrations should be determined by the end user.

性 状 Liquid

浓 度 2mg/1ml

免 疫 原 Native Goat IgG

亚 型 IgG

纯化方法 affinity purified by Protein A

缓 冲 液 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

注意事项 This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品介绍 Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

### 【备注】

本产品仅供科研使用。在确认产品质量出现问题时，本公司承诺为客户免费更换等量的质量合格产品。