

## M5 Donkey Anti-rabbit IgGFITC antibody

Product	Unit	Cat.#
M5 Donkey Anti-rabbit IgGFITC antibody	100ul	MF757-01

【Storage】: Store at -20℃ for one year. Avoid repeated freeze/thaw cycles.

### 【产品介绍】:

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

### 【产品参数】:

抗体来源 Donkey

克隆类型 Polyclonal

交叉反应 Rabbit

产品应用 Flow-Cyt=1:100-1000 IF=1:100-1000

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

性 状 Liquid

浓 度 2.0 mg/ml

免 疫 原 Native rabbit IgG

亚 型 IgG

纯化方法 affinity purified by Protein G

缓 冲 液 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

**Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"**