

M5 Goat Anti-Rabbit IgG/Alexa Fluor 594 antibody 使用说明书

Product	Unit	Cat.#
M5 Goat Anti-Rabbit IgG/Alexa Fluor 594 antibody	100ul	MF937-01

保存条件: Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

中文名称: AF594 标记的山羊抗兔 IgG H&L

别名: Goat Anti-Rabbit IgG H&L (AF594); Immunoglobulin G;

抗体来源: Goat

克隆类型: Polyclonal

交叉反应: Rabbit

产品应用: IF=1:200-1000, Flow-Cyt=1:50-200
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

性状: Liquid

浓度: 2.0 mg/ml

免疫原: Native rabbit IgG

亚型: IgG

纯化方法: affinity purified by Protein G, nonspecific adsorbed

缓冲液: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

注意事项: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.