

## M5 HiPer 重组人 VEGF165 使用说明书

产品名称	单位	货号
M5 HiPer 重组人 VEGF165	10µg	MFP06008-01
M5 HiPer 重组人 VEGF165	50µg	MFP06008-05

### 【储存条件】

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks.

Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

### 【质量控制】

Greater than 95% as determined by reducing SDS-PAGE.

### 【产品浓度】

100 µg/ml.

### 【产品简介】

Human Vascular endothelial growth factor (VEGF), also known as VEGF-A and vascular permeability factor (VPF), belongs to the platelet-derived growth factor family of cysteine-knot growth factors. It is a potent activator in vasculogenesis and angiogenesis both physiologically and pathologically. VEGF-A has 8 differently spliced isoforms, of which VEGF165 is the most abundant one. VEGF165 is a disulfide-linked homodimer consisting of two glycosylated 165 amino acid polypeptide chains. VEGF stimulates the cellular response through binding to tyrosine kinase receptors VEGFR1 and VEGFR2 on the cell surface. It is widely accepted that VEGFR2 mediate almost all of the known cellular responses to VEGF while the function of VEGFR1 is less defined and is thought to modulate the VEGFR2 signaling.

### 【注意事项】

Always centrifuge tubes before opening. Do not mix by vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100 µg/ml.

Dissolve the lyophilized protein in ddH<sub>2</sub>O.

Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

### 【备注】

本产品仅供科研使用。在确认产品质量出现问题时，本公司承诺为客户免费更换等量的质量合格产品。