




## Product Specification

Product Name: **Agarose**

Product Number: **MF103**  
CAS Number: **9012-36-6**



### TEST

### Specification

|                                    |                           |
|------------------------------------|---------------------------|
| Solubility (Colour)                | Colorless to Light Yellow |
| Solubility (Turbidity)             | Clear to Hazy             |
| 1.5 g plus 100 mL H <sub>2</sub> O |                           |
| Congearing Temperature             | 34.5 - 37.5 °C            |
| Loss on Drying                     | ≤ 10 %                    |
| Appearance (Colour)                | White to Off-White        |
| Appearance (Form)                  | Powder                    |

Specification: PRE.0.ZQ5.10000002515

Mei5bio warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [mei5bio.com](http://mei5bio.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.